



## White / Black RF Face Lifting Machine 5 Modes Radio Frequency Skin Tightening Equipment

### Our Product Introduction

For more products please visit us on [hifubeauty-machine.com](http://hifubeauty-machine.com)

#### Basic Information

- Place of Origin: China Guangzhou
- Brand Name: GLM
- Certification: CE
- Model Number: X10
- Minimum Order Quantity: 1
- Price: confer
- Packaging Details: Standard export packing
- Delivery Time: 5-8days
- Payment Terms: L/C, D/A, D/P, T/T, Western Union, MoneyGram
- Supply Ability: 100 Piece/Pieces per Month



#### Product Specification

- Name: RF Face Lifting Machine
- Number Of Modes: 5
- Feature: Skin Whitening
- Led Wavelength: 620nm, 415nm, 530nm
- Material: ABS + Stainless Steel
- Rf Frequency: 1MHz
- Output Power: 10W
- Other Name: Radio Frequency Facial Machine
- Color: White, Black
- Name Lifting Body Shaping Machine: 360 VELA-SHAPE Slimming Machine, Commercial VELA-SHAPE Slimming Machine, White RF Face Lifting Machine
- Application For Commercial Warranty 2 Years OEM&ODM



#### More Images





## Product Description

### Product Description:

This Fiber Optic Instrument features multiple LED wavelengths, including 620nm, 415nm, and 530nm, which work together to help reduce the appearance of fine lines and wrinkles, while also brightening and rejuvenating the skin.

One of the key benefits of this EMS Beauty Instrument is that it can help with skin whitening, making it an excellent choice for those looking to improve their skin tone and achieve a more even complexion. In addition, it is also effective for face slimming, face lifting, and anti-aging.

Available in both white and black, this CET RET RF Machine is easy to use and can be incorporated into any skincare routine. Whether you are looking to address specific skin concerns or simply want to maintain a healthy, youthful glow, this device is a great choice.

### Features:

Product Name: RF EMS Beauty Instrument

Output Power: 10W

Other Name: Radio Frequency Facial Machine

Function: Face Slimming, Face Lifting, Anti-aging

Number Of Modes: 5

Led Wavelength: 620nm, 415nm, 530nm

This Face Beauty Instrument, also known as the Radio Frequency Facial Machine or Elight Rf Laser Machine, has an output power of 10W and 5 different modes to choose from. It is designed for face slimming, face lifting, and anti-aging treatments. Additionally, it features LED lights with wavelengths of 620nm, 415nm, and 530nm. This Fiber Optic Instrument is a must-have for anyone looking to improve their skin's appearance.

### Technical Parameters:

Output Power	10W
Number Of Modes	5
Led Wavelength	620nm, 415nm, 530nm
Function	Face Slimming, Face Lifting, Anti-aging
Rf Frequency	1MHz
Voltage	100-240VAC
Material	ABS + Stainless Steel
Color	White, Black
Other Name	Radio Frequency Facial Machine
Treatment Area	Face, Neck, Body

This is a technical table for a Fractional RF Beauty Machine, a type of Beauty Care Instrument or Beauty Medical Equipment.





handle 1 for face



handle 2 for body



handle 3 for body



handle 4 for body

### Applications:

The RF EMS Beauty Instrument operates at a frequency of 1MHz, making it an efficient and effective tool for skin rejuvenation. This product is perfect for people who are looking for an OEM Thermage RF Machine. It works by emitting radiofrequency waves that penetrate deep into the skin, stimulating collagen and elastin production, which in turn tightens and lifts the skin.

The RF EMS Beauty Instrument is a versatile product that can be used on the face, neck, and body. It is ideal for people who are looking for Skin Tightening Slimming Beauty Equipment. Its EMS function helps improve blood circulation, which can reduce the appearance of cellulite and aid in slimming. This product is also perfect for those who would like to have a VELA-SHAPE Slimming Machine.

One of the key features of the RF EMS Beauty Instrument is its ability to whiten the skin. It is perfect for people who want to achieve a brighter and more even skin tone. This product is perfect for people who want to have a radiant glow. The RF EMS Beauty Instrument has a voltage range of 100-240VAC, making it suitable for use in different countries. This product is perfect for salon owners, beauty clinics, and individuals who want to have their own beauty equipment. The RF EMS Beauty Instrument is perfect for various occasions and scenarios. It is perfect for people who want to look their best for special events such as weddings, parties, and other social gatherings. It is also perfect for people who want to maintain their youthful appearance. This product is also perfect for people who want to relax and pamper themselves after a long day at work or during their free time.

### Customization:

Customize your RF EMS Beauty Instrument with our product customization services!

Originating from Guangzhou, China, our Radio Frequency Facial Machine offers 5 different modes and an output power of 10W. It operates on a voltage range of 100-240VAC and offers three different LED wavelengths: 620nm, 415nm, and 530nm.

Along with these standard features, we also offer product customization services to make your Fiber Optic Instrument truly unique.

Whether you want to add extra LED wavelengths, adjust the output power, or modify the device's appearance, we can help.

Our services are perfect for beauty salons and spas looking to offer the latest in skin rejuvenation technology. With our customization options, you can create your own Salon Thermage RF Machine or CET RET RF Machine that perfectly meets the needs of your business and customers.



**Guangzhou Lizhirun Technology Development Co., Ltd**



18688905627



cenglinlin25@gmail.com



hifubeauty-machine.com

Room 502, 5th Floor, Building F, Huixin Industrial Zone, No. 89, Hetai Road, Hebain Village, Baiyun District,  
Guangzhou